

What are the benefits?

- No risk of flooding from running water
- Eliminate the cost of water purchase and disposal
- For solvent volumes from 5 ml up to 1 litre
- · Helps meet sustainable water reduction targets

Average water consumption of common household appliances compared to a Findenser and water condenser



Findenser[™]- prevents flooding and saves water

Replaces water-cooled condensers in over 95% of common chemistry applications

How does Findenser work?

- Findenser comprises an internal glass condenser and an external, finned aluminium jacket, between which a small amount of water is permanently sealed.
- The glass condenser design has a greater internal surface area than traditional air condensers, increasing heat transfer capacity.
- The finned jacket fits around the glass condenser, further increasing the external surface area.
- The result is a 'SUPER air condenser'.



Findenser requires no running water to operate. Water is a precious resource. It makes little economic or environmental sense to waste thousands of litres just to cool a single condenser.

Performance testing

A range of solvents, in identical flasks and set-ups, were tested with a Findenser, water condenser and air condenser to record solvent loss by weight.

Findenser compared to an air condenser

For synthesis with low boiling point solvents, Findenser showed a significant improvement in solvent retention. With acetone or DCM the reaction boiled dry when using an air condenser, vet Findenser retained 95% of the solvent under the same conditions.

For synthesis with medium boiling point solvents, Findenser delivered improved solvent retention particularly with larger volumes and high temperatures.

Findenser compared with a water condenser

Under identical conditions, a standard Findenser retained solvent to the same level as a water condenser (with the exception of diethyl ether).



Radleys Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom

Choice of B14, B19, B24 and B29 joint sizes







B19